Indications:

Absorbest BioShield is a sterile superabsorbent wound dressing for use on exuding wounds.

Product properties:

Absorbest BioShield has a superabsorbent core which absorbs and retains exudate. The exudate is absorbed vertically by the dressing from the white wound contact side and is retained in the core forming a gel, thus minimizing the risk of maceration to the peri-wound skin while maintaining a moist wound environment. The white wound contact layer can have a preventive effect on biofilm formation on its surface and immediate surroundings, as shown *in vitro*. When exudate is absorbed, the wound contact layer surface becomes hydrogel-like and moist. Absorbest BioShield has a hydrophobic green backing layer (side up) to minimize moisture strike-through.

Instructions for use:

Absorbest BioShield should be used on exuding wounds (patients from 1 year and older) and handled by healthcare professionals. The dressing can be used as either primary or as secondary dressing.

A) Clean the wound as per local protocol. Dry the surrounding skin carefully. **B)** Choose a size of dressing depending on the wound size; the core edge of the dressing should overlap the wound with a 2 cm margin as the sealing edges do not absorb exudate. Apply the dressing with the white side in direct contact with the wound bed or the primary dressing. **C)** Secure the product in place. The dressing can be used in combination with compression therapy. **D)** The dressing should be changed according to the clinical condition of the wound, or before saturation. **E)** Dressing removal can be facilitated by moistening the dressing. **F)** Dispose of the dressing as clinical waste according to local regulations.

If the wound deteriorates, consult relevant healthcare professionals. In case a serious adverse event has occurred in relation to usage of the device, this should be reported to Absorbest AB and the competent authority in the country in question.

Caution:

Do not cut the dressing, the cover must remain intact. Do not use the dressing if the sterile packaging is damaged. Single use product, do not re-sterilize, can not be re-used. Important to use the white wound contact side towards the wound! Absorbest BioShield should not be used on dry wounds.

Contraindications:

Absorbest BioShield dressing should not be used on eyes, mucous membranes, tendons or in wound cavities because the dressing swells during absorption. Absorbest BioShield dressing should not be used on persons with a known sensitivity to the dressing.

-  Do not reuse. Reusing will cause cross contamination.
-  Do not use if the sterile barrier or packaging is compromised.
-  Do not cut.

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